

PREVENTIS *SmarTest Pro*® Vitamin D

Rapid smart test for
the determination of
vitamin D level



A modern solution: The combination of a rapid test and a smartphone app for the accurate determination of vitamin D levels in minutes.

The rapid test uses a small blood sample from the fingertip and the test result is displayed directly on the smartphone. Performing the test in the doctor's office or at a pharmacy is simple and makes immediate medical counseling possible.

Preventis is your renowned partner for health care. Our research provides you today this modern smartphone-driven rapid test for the reliable determination of vitamin D.

Vitamin D smart test



Product Code	Description
PST80100G5	SmarTest Pro Vitamin D smart test for 25(OH)Vitamin D in finger-prick blood samples. 5 tests.


BIOHIT HealthCare Ltd
Pioneer House, Pioneer Business Park
North Road, Ellesmere Port
CH65 1AD, United Kingdom

Tel: +44 151 550 4 550
sales@biohithealthcare.co.uk

www.biohithealthcare.co.uk


BIOHIT HealthCare
Innovating for Health

BIOHIT HealthCare
Innovating for Health



PST80100G5EN2019/04

Vitamin D: vital for health



Vitamin D plays an important role in our metabolism, especially for bone health. The body can produce vitamin D in the skin with the help of sunlight. Therefore, it is often called the "sunshine vitamin". Unfortunately, reduced sunlight hours in the winter, lack of time outdoors, and use of sunscreen can easily lead to a vitamin D deficiency.

Ideal vitamin D level

- Promotes bone health and prevents osteoporosis or rickets
- Protects against infections, for example cold


Vitamin D deficiency and its relation to cardiovascular disease, autoimmune diseases, and even cancer and diabetes are also discussed in medical science.

Tips for achieving an optimal level of vitamin D

- ✓ Frequent outdoor activities (15–20 minutes of sun with more than 40 % of skin exposure is optimal)
- ✓ Eating vitamin D rich foods: fish (salmon, herring), liver, cheese, mushrooms and avocados
- ✓ Adding vitamin supplements to diet

PREVENTIS *SmarTest Pro*® Vitamin D

How does SmarTest Pro® work?



The test is performed by professionals, for example, in the doctor's office or at a pharmacy:

- ✓ Painlessly collect a blood sample from the fingertip
- ✓ Apply the sample to the test device
- ✓ Incubate test for 15 minutes
- ✓ Scan the test device with the smartphone camera app*
- ✓ The **Preventis SmarTest Pro® Vitamin D** app calculates the vitamin D level
- ✓ The result is displayed immediately on the smartphone of the doctor or the pharmacist

For more information visit:

www.biohithealthcare.co.uk/smartestpro-vitd

PREVENTIS *SmarTest Pro*® Vitamin D

What are the advantages of SmarTest Pro®?

Quickly detect vitamin D deficiency and continuously monitor optimal levels

Quantitative Test results are available within a few minutes

- ✓ Performing the test is simple, for example in the doctor's office or at a pharmacy
- ✓ Painless collection of a blood sample from the fingertip
- ✓ Medical counseling available immediately following the test
- ✓ Optimal adjustment of vitamin D supplementation based on continuous monitoring

When taking vitamin D supplements it is advisable to check the effectiveness on a regular basis.

Preventis SmarTest Pro® Vitamin D app required

